

Wellbeing: Anti-Ligature Grab Rail - Stainless Steel - 965 x 600mm Toilet Assisted 90 ° Left Hand

These stainless steel grab rails are ,anti-ligature,, which means that any points where a cord, rope or bed sheet can be looped or tied through to create a point of ligature are eliminated.

These anti-ligature grab rails are ideal for use in hospitals, correctional centres and mental health facilities where self-harm is a concern.

KEY FEATURES

- Fully welded construction in grade 304 38mm stainless steel tubing.



PRODUCT CODE

| | |
|----------------|---|
| ECALIG965L90SS | Wellbeing: Anti-Ligature Grab Rail - Stainless Steel - 965 x 600mm Toilet Assisted 90 ° Left Hand |
|----------------|---|

Due to ongoing Research and Development, specifications may change without notice.

Component specifications may change on some export models.

Refer to warranty statement for warranty details - www.enware.com.au/warranty.

Fixing screws are not supplied with product.

Ensure grab rail is installed on a solid wall surface with firm wall backing using appropriate fixing method.

Structural engineers or licensed tradespersons should be consulted for appropriate advice on suitability of the fasteners and wall structure before affixing or installing any Enware products.

Enware does not provide any warranty in relation to bolts, screws or other fasteners used with Enware products.

TECHNICAL DATA

| | |
|-------------|---|
| Load Rating | Enware Grab Rails that are manufactured within the guidelines of AS1428 are compliant with being able to withstand 1100N of force applied in any direction. |
|-------------|---|

